

CCDC153-coiled-coil domain containing 153 Gene

Reference: PTXBC101443

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	coiled-coil domain-containing protein 153; coiled-coil domain containing 153		
For research use only.			

GENE DESCRIPTION	coiled-coil domain containing 153
GENE ID	283152
ORIGIN SPECIES	Human
ACCESSIONS	BC101443
VECTOR	pDONR223
NCBI SYMBOL	CCDC153
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtcgccagtgcctgcaggagatgcaaacccacagcaagcagctggaggaagaagtcaaggcctcggggcagctagagcatgccaaaggaggctgcagctgccgggaagaggctgaacaagctctggagagcgggaccaggccctggctcagctt
cgggccacatggcagacatggagcgaagtatgaggaaatcttacgacagctctggacaggctctggccaagctgagagccatcaagcagcagtgaggatggggcagcactgagacttcacgccaggcacaaggagcagcaagccagtttgactacccccctggatcttg
aggccaccagccctagctctag