

SSX3-synovial sarcoma, X breakpoint 3 Gene

Reference: PTXBC103862

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Kruppel-associated box containing gene product SSX3; protein SSX3; CT5.3; cancer/testis antigen 5.3; synovial sarcoma, X breakpoint 3; SSX family member 3		

For research use only.

GENE DESCRIPTION	synovial sarcoma, X breakpoint 3
GENE ID	10214
ORIGIN SPECIES	Human
ACCESSIONS	BC103862
VECTOR	pDONR223
NCBI SYMBOL	SSX3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacggagatgacaccttgcaggagaccacggttggtcctcaataaccagagagatacaaaagcccttcgatgatattgccaaatactctctaaggaagagtgaggaaaagatgaaagtctcggagaaaatcgtctatgtgtatatgaagaaaagtatgaggccatgactaaact
agggttcaaggccatcctccatcttcatgctgaataaacgggtcacagactccagggggaatgatttgataatgaccctaaccgtgggaatcaggttcaacgtctcagatgactttggcaggctccagggaatctccgaagatcatcccaagaagccagcagaggaaggaatgt
ttcgaagggaagtgcagagcatctgcccacaaaacgatgggaacagctgtgcccccgggaaaacaaactacctctgagaagattaacatgatatctggacccaaaaggggggaacatgcctggaccacagactgcgtgagagaaagcagctggtgatattgaagagatcagc
gatcctgagggaatgatgagtaa