



Reference: PTXBC103911

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TCEB3B; HsT832; TCEB3L; RNA polymerase II transcription factor SIII subunit A2; transcription elongation factor (SIII) elongin A2 elongin A2; transcription elongation factor B polypeptide 3B; transcription elongation factor B subunit 3B; elongin A2		

For research use only.

GENE DESCRIPTION	transcription elongation factor B polypeptide 3B (elongin A2)
GENE ID	51224
ORIGIN SPECIES	Human
ACCESSIONS	BC103911
VECTOR	pDONR223
NCBI SYMBOL	TCEB3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]