



Reference: PTXBC103926

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PCDH-GAMMA-B1; protocadherin gamma-B1; protocadherin gamma subfamily B, 1		
For research use only.			

GENE DESCRIPTION	protocadherin gamma subfamily B, 1
GENE ID	56104
ORIGIN SPECIES	Human
ACCESSIONS	BC103926
VECTOR	pDONR223
NCBI SYMBOL	PCDHGB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr