

KRTAP13-3-keratin associated protein 13-3 Gene

Reference: PTXBC103965

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP13.3; keratin-associated protein 13-3; keratin associated protein 13-3		
For research use only.			

GENE DESCRIPTION	keratin associated protein 13-3
GENE ID	337960
ORIGIN SPECIES	Human
ACCESSIONS	BC103965
VECTOR	pDONR223
NCBI SYMBOL	KRTAP13-3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctacaactgtgtcttagaaactctctctgtctccacgggggttacttgcaactaccaggtctctctgtgtggtctctctacccagcaacctgggtctacagcaactgtctctccagcacctgccagctgggtctctctctatagggctgtcaggagacctgtggagg
cccaacagctgtcagacattgtgtgtgagtcagccctgccacacctctgtactacccaggaactcacatgctctgcaattctgtcctgactatgcatgttgggtctcgggttttgatccaatagctgtctccctgagctgtggatccaggagctgtctctcactggctgtggatcc
aatggcttcagatatctgaattatagaatccatacctccccctccagagttatagatccagattctgcatccaatctatttccacctagaagggtgttcattcatctgttatcagccattctgtatctgttttctactga