



Reference: PTXBC105640

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PCDH-BETA10; PCHB10; protocadherin beta-10; PCDH-beta-10; protocadherin beta 10		

For research use only.

GENE DESCRIPTION	protocadherin beta 9
GENE ID	56127
ORIGIN SPECIES	Human
ACCESSIONS	BC105640
VECTOR	pDONR223
NCBI SYMBOL	PCDHB9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]