

## ZNF182-zinc finger protein 182 Gene

Reference: PTXBC106875

|                |                                                                                                                   |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                               | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                              | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | HHZ150; KOX14; ZNF21; Zfp182; zinc finger protein 182; zinc finger protein 21 (KOX 14); zinc finger protein KOX14 |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | zinc finger protein 182                                                                          |
| <b>GENE ID</b>          | 7569                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC106875                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | ZNF182                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgacaccagcatctgcatctggtgagactcaggaagcttctactatggcagaaggcaagagggaacaggggctagtgcatttgaagatgtagctgtgatttcacccaggaggagtgccagctgaacccaccacagaggaccctgtacagagacgtgatctggagaccta  
cagcaacttggtcttggggcagcaggttaccaaaccaaacctcatctcaagttggaggtagaagaatgccggcagaaggaaaaatccatttggaaactttccagaagctctgcaagttgatgaacagattgagaggcaacatcaggatgaccaagataaatgtctgtatgcaagt  
tggatttttgacaagaaaacaattatccaagagtgctcgtgactgtcatgatttggaaacatacttcatctgagtacaaccttggcttcaatatacaagaccgataaacacgaatcatttggaaataatggttagataatttagacttatttagagaagtctgcagaaaaataatgat  
aatggatgtgcacaaatgttctccatactgagtatgagaaaaacaatcctggaaatgaagccctatggctataaagagtggtgggaaaggtcttaggcgaagaaagggccttagtctacatcagagaattaaaaatggagagaaccccttgaatgtactcatgtaggaaaaccttcagcaaga  
agtacacctcatgtacattggagaactcagggagagaaaccccttggatgtaccgaatgtggaaaagctttagccaaaaatctcagctcattatacacttgcgaactcatacaggagagagaccctttagtgcctgaatgtggaaaagcccttcagagaaaagtcaactgtcatcaca  
attacagactcacacagagagaaaaacccctatgaatgaatgtgaaaagcccttcactcagaagttccaacctcattgtccatcagaaaacccacactggagagaaaacccctatgaatgcataaatgtggagaatctttacacagaagcttgatctaattatacatagaccataca  
ggaaagaaaccccatgaatgaatgagtgtaagaaaacttfcagtgataagtcacactcattatacaccagagaactcatacggagagagaacccataaatgtactgtgaatgtgggaagcttctcaatgagaagtcacccctcattgtgcataaagaactcatacaggagagaacccctatg  
aatgtgatgtgtggaaaaaccttcagcaaaagtcaaaccttggctacatcagagaactcattcaggagagaaccccttgaatgaatgaatgtggaagcccttctcagaagtcctacatcgtgcacatcagagaggtcatacaggagagaacccctacgagtgcaatgaatgtga  
aaaagcattttcagagaatcattatcattatacacaagaactcatacagaagaaaaacccctataaatgaatgaatgtggcaagcccttcagagaaaaagtcacaaagctcattatacagagaattcatacaggagagaacccctatgaatgtcctgtgtgtggaagccttttagccagaa  
gtcacagctcataatacatcagagaacgcacacgggagagaacccctatgcctgactgagtggtggcaagcccttcgagaaaaagtcacattcgtacatcaagaactcactggagagaacccctataaatgtacagaatgtgggaagcccttaccacaaaatcaaaccttattgt  
acatcagcgaacccatcagagaaagaaagcccatggaagagccacactcgaagtcgaagtcacacattg