

HIST2H2BA-histone cluster 2, H2ba Gene

Reference: PTXBC106916

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	histone cluster 2, H2ba
GENE ID	337875
ORIGIN SPECIES	Human
ACCESSIONS	BC106916
VECTOR	pDONR223
NCBI SYMBOL	HIST2H2BA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcatcatgaactccttcgtcaacgacatctcgagcgcatcgccggagaggcgtcccgcctggcgcaactacaacaagcgtcccatcacatccgtgagatccagacggccgtgcgcctgctgctgcccgagctggccaagcacgccgtgtccgagggcaccagggc
ggtcaccaagtacaccagctcgaacattttctag