



Reference: PTXBC106937

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 791		
For research use only.			

GENE DESCRIPTION	zinc finger protein 791
GENE ID	163049
ORIGIN SPECIES	Human
ACCESSIONS	BC106937
VECTOR	pDONR223
NCBI SYMBOL	ZNF791
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]