

C11orf85-chromosome 11 open reading frame 85 Gene

Reference: PTXBC106952

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C11orf85; membrane-anchored junction protein; membrane anchored junction protein		
For research use only.			

GENE DESCRIPTION	chromosome 11 open reading frame 85
GENE ID	283129
ORIGIN SPECIES	Human
ACCESSIONS	BC106952
VECTOR	pDONR223
NCBI SYMBOL	C11orf85
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttttaaaccctttacccacccgtttccagagacgaggtttcttcacgagcccaatgtgtataaattcaaatcagatatgggaagagtatcagaggagaagagatagaaataaggagtcacccaggagctggaggattctgtccgcgtggtcttgggaaactggacaatc
ttcagcccttgcacagaaacttcattgtattccctataaaagcaaatgggagagagtttccacctgaaattcaaacatggggaattatctgtatccctacccattgttttactctatatgtggagatgaaatggttcattgaaacctgtccctgggaaaccaataagtgcagtcctctt
gggttgggtccagttgagaagaagcagtaggagctgtgatgaggaaacgaaacacatggacgagccagctcccccagcaggccagggctggacagaataggaaagaaaaacccaaggattgcaggagactctggcctctgatactactgtccagaacaagattctg
agtggggacacagcctgccaggcgaattgtccacccctgcagcacaactcacctccactaaggagcgagcagccacgggtcttcttgggtttctaa