

PCGF3-polycomb group ring finger 3 Gene

Reference: PTXBC107061

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DONG1; RNF3; RNF3A; polycomb group RING finger protein 3; RING finger protein 3A; ring finger protein 3; polycomb group ring finger 3		

For research use only.

GENE DESCRIPTION	polycomb group ring finger 3
GENE ID	10336
ORIGIN SPECIES	Human
ACCESSIONS	BC107061
VECTOR	pDONR223
NCBI SYMBOL	PCGF3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgaccaggaagatcaagctgtgggacatcaacgccacatcacctgccgctgtgcagcgggtacctcatcgacgccaccacgtgaccgagtgctgtgcacacctctgcaggagctgcctgggtgaagtaacctggaggagaacaacacctgccccacctgcaggattgtgateca
ccagagccacccctgcagtagatcgggtcatgacagaacctgcaagatattgtttacaaattgtaccaggcctcagaagcggaaatgagaaagcagagggtgtctatcacaattggcctgaggtgcccggagacatcaagggggagacctgctctgcaaaacagcacttag
attccatcggaatggtgaaacaaagcagacgacagtcaaaacaaagaggccgaggagagagcggagagacaacgactaccaccgcagcgcagcagcaggtgagcatctgcctggagtgtacagcagcaaacctgcgcgggctgaagcgggaagtggatccgctgctcag
cccaggcgaccgtctgcatctgaagaagttcatgcgcaaaaaactcaacctttcatcttaacgagctggacattttatgcaacgaggagatcctggcaggaccacacactcaagttcgtggtgtcactaggtggagattcaagaaggcggcctctgctgactacagaccaaga
tggaactgctgtga