

TMEM95-transmembrane protein 95 Gene

Reference: PTXBC107110

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UNQ9390; transmembrane protein 95		
For research use only.			

GENE DESCRIPTION	transmembrane protein 95
GENE ID	339168
ORIGIN SPECIES	Human
ACCESSIONS	BC107110
VECTOR	pDONR223
NCBI SYMBOL	TMEM95
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggaggctggcactaggcgggttttctggcagccgccaggctgtgtctctgtcgcctccagccacgaactgtcaggccgcctggctcgtcgcagccagatggaggccaggcagaagggaatgtggggcctcccgactctcggcctttgccttagatgaggtgtcca
tgaacaaagtacagagaagactcacagagtcctgagggtcatggaatcaaaaggctgtctcctcactccctcatattggagttggttcgaaagaccaagctccctgagtagaccagggaagctctctgtccccgcctgccggggcagcaccacgctgtacaactgtccacctg
caaggggacggaggtgtctgtctgccccgaaagcgtgtctccagggaagtcaggatctttgggaagccaagattctgtcctcctcctctcggagcttctcgtctctggtgtgtctgagcctcctggtggagtcaccacacccaagcaaaagtggcttga