

C6orf184-chromosome 6 open reading frame 184 Gene

Reference: PTXBC107113

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	chromosome 6 open reading frame 184
GENE ID	221261
ORIGIN SPECIES	Human
ACCESSIONS	BC107113
VECTOR	pDONR223
NCBI SYMBOL	C6orf184
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcaaaacttgcagaaaaggcagcttttaaaccttaagcgaaactctacaactgatatgctctctgcatgatctgtgacctacgttgcagttttgctaagcttggcaattgtccagcatgtttgaatttctcgaagccaacctttgagaggtgactggggagggaactgagggcat
tgggtctgagcttcaagagctgcagaacatgattgacagctccagagcccccaagacctatccgggtggccaggcactctctccggaggaggttatattttgcagtttgacgctgcagtaaggcatctcatccgaagaacattttggcagctggaaatgttctgctaccagctct
gtcagagacggcatgtgcatgggtctaccagcactgagcaactctctcaggaagacattttgcttcacagctcagctgccccagccactggatccacggagcctccaggcatttgagctgtttccttggagagcatttctggaagatggaggaccattccagttatgagtaacagccc
agataccctagaatataatgtaggttaggtag