



Reference: PTXBC107168

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C1orf112; chromosome 1 open reading frame 112		
For research use only.			

GENE DESCRIPTION	chromosome 1 open reading frame 112
GENE ID	55732
ORIGIN SPECIES	Human
ACCESSIONS	BC107168
VECTOR	pDONR223
NCBI SYMBOL	C1orf112
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]