



Reference: PTXBC109125

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CLR19.4; HYDM; NALP7; NOD12; PAN7; PYPAF3; NACHT, LRR and PYD domains-containing protein 7; NACHT, LRR and PYD containing protein 7; NACHT, leucine rich repeat and PYD containing 7; PYRIN-containing Apaf1-like protein 3; nucleotide-binding oligomerization		

GENE DESCRIPTION	NLR family, pyrin domain containing 7
GENE ID	199713
ORIGIN SPECIES	Human
ACCESSIONS	BC109125
VECTOR	pDONR223
NCBI SYMBOL	NLRP7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1