



Antibodies and Proteins

## ProteoGenix SAS

## Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

[www.ProteoGenix.science](http://www.ProteoGenix.science)

## CCDC140-coiled-coil domain containing 140 Gene

Reference: PTXBC109388

|                               |                                                                              |                     |                |
|-------------------------------|------------------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                                          | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                                         | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | coiled-coil domain-containing protein 140; coiled-coil domain containing 140 |                     |                |
| <b>For research use only.</b> |                                                                              |                     |                |

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | coiled-coil domain containing 140                                                                 |
| <b>GENE ID</b>          | 151278                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC109388                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CCDC140                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

[illegible]

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)