

DEFB127-defensin, beta 127 Gene

Reference: PTXBC110085

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf73; DEF-27; DEFB-27; DEFB27; bA530N10.2; hBD-27; beta-defensin 127; beta-defensin 27; defensin beta 127		

For research use only.

GENE DESCRIPTION	defensin, beta 127
GENE ID	140850
ORIGIN SPECIES	Human
ACCESSIONS	BC110085
VECTOR	pDONR223
NCBI SYMBOL	DEFB127
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggctctcatgacattgcaattctgctgtccagaaacccacagtaaccgaacaacttaagaagtgctggaataactatgtacaagacattgcaggaaaatctgcagagtaaatgaagtcctgaggcactatgtgaaatgggagatctgtgcctcaatatcaaggaaactggaag
catgtaaaaaaattacaagccacctagtccaaagccagcaacacttgcactgactctcaagactatgttacaataatagaaaatttccaagcctgaagacacagctctacataa