



Reference: PTXBC110542

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HZF3; zinc finger protein 268; zinc finger protein 3; zinc finger protein HZF3		
For research use only.			

GENE DESCRIPTION	zinc finger protein 268
GENE ID	10795
ORIGIN SPECIES	Human
ACCESSIONS	BC110542
VECTOR	pDONR223
NCBI SYMBOL	ZNF268
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr