



Reference: PTXBC111413

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SCDO2; bHLHc6; mesoderm posterior protein 2; class C basic helix-loop-helix protein 6; mesoderm posterior 2 homolog; mesoderm posterior basic helix-loop-helix transcription factor 2; mesoderm posterior bHLH transcription factor 2		

For research use only.

GENE DESCRIPTION	mesoderm posterior 2 homolog (mouse)
GENE ID	145873
ORIGIN SPECIES	Human
ACCESSIONS	BC111413
VECTOR	pDONR223
NCBI SYMBOL	MESP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]