



Reference: PTXBC112199

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ERBB4 transcript variant I20DEL; ERBB4 transcript variant I12DEL; ALS19; HER4; p180erbB4; receptor tyrosine-protein kinase erbB-4; avian erythroblastic leukemia viral (v-erb-b2) oncogene homolog 4; human epidermal growth factor receptor 4; proto-oncogene-		

For research use only.

GENE DESCRIPTION	v-erb-a erythroblastic leukemia viral oncogene homolog 4 (avian)
GENE ID	2066
ORIGIN SPECIES	Human
ACCESSIONS	BC112199
VECTOR	pDONR223
NCBI SYMBOL	ERBB4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]