



Reference: PTXBC112238

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NSMASE2; sphingomyelin phosphodiesterase 3; nSMase-2; neutral sphingomyelinase 2; neutral sphingomyelinase II		

For research use only.

GENE DESCRIPTION	sphingomyelin phosphodiesterase 3, neutral membrane (neutral sphingomyelinase II)
GENE ID	55512
ORIGIN SPECIES	Human
ACCESSIONS	BC112238
VECTOR	pDONR223
NCBI SYMBOL	SMPD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr