

HIST1H2AI-histone cluster 1, H2ai Gene

Reference: PTXBC112254

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2A/c; H2AFC; histone H2A type 1; H2A histone family, member C; H2A.1; histone 1, H2ai; histone cluster 1, H2ai; histone cluster 1 H2A family member i		

For research use only.

GENE DESCRIPTION	histone cluster 1, H2ai
GENE ID	8329
ORIGIN SPECIES	Human
ACCESSIONS	BC112254
VECTOR	pDONR223
NCBI SYMBOL	HIST1H2AI
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggcggtggcaagcaggaggcaaaagctcgcgcaaggccaagaccgctcttctcgggcccggcttcagttcccgtaggcccagtgcatcgctgctccgcaaaggcaactatgcggagcgggtcggtgctggagcgcgggtgtacctggcggcggtgctggagctacgtg
accgccgagatcctggagctggctggcaacgcggcgccgacacaagaagactcgcacatcccgctcactcagctggccatccgcaacgatgagagctcaacaagcttctgggcaaagtcaccatgcacagggtggcgctcgtgcccaacatccaggccgtgctactgccc
aagaagaccgagagccaccacaaggcgaagggaagtag