



Reference: PTXBC112307

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPR133; PGR25; adhesion G-protein coupled receptor D1; G-protein coupled receptor 133; G-protein coupled receptor PGR25; adhesion G protein-coupled receptor D1		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 133
GENE ID	283383
ORIGIN SPECIES	Human
ACCESSIONS	BC112307
VECTOR	pDONR223
NCBI SYMBOL	GPR133
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]