

LALBA-lactalbumin, alpha- Gene

Reference: PTXBC112316

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LYZG; alpha-lactalbumin; lactose synthase B protein; lysozyme G; lysozyme-like protein 7; lactalbumin alpha		

For research use only.

GENE DESCRIPTION	lactalbumin, alpha-
GENE ID	3906
ORIGIN SPECIES	Human
ACCESSIONS	BC112316
VECTOR	pDONR223
NCBI SYMBOL	LALBA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggtctttgtccctctgttctgtgtggcctcctgttccctgcctcctgccaagcaattcacaaatgtgagctgtccagctgctgaaagacatagatggttatggagcctccttgcctgaattgatctgtaccatgtttcacaccagtgggttatgacacacaagccatagttgaaa
caatgaaagcacggaatatgactcttccagatcagtaataagctttgtgcaagagcagccaggtccctcagtcaggacatctgtgacatcctgtgacaagttcctgcatgatgacattactgatgacataatgtgtgccaagaagatcctggatattaaagggaattgactactggttg
ccataaagccctctgcactgagaagctggaacagtggctttgtgagaagttgtga