



Reference: PTXBC113372

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PCDH-ALPHA12; protocadherin alpha-12; KIAA0345-like 2; PCDH-alpha-12; protocadherin alpha 12		

For research use only.

GENE DESCRIPTION	protocadherin alpha 12
GENE ID	56137
ORIGIN SPECIES	Human
ACCESSIONS	BC113372
VECTOR	pDONR223
NCBI SYMBOL	PCDHA12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]