



Reference: PTXBC113407

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NUCLING; uveal autoantigen with coiled-coil domains and ankyrin repeats; nuclear membrane binding protein		

For research use only.

GENE DESCRIPTION	uveal autoantigen with coiled-coil domains and ankyrin repeats
GENE ID	55075
ORIGIN SPECIES	Human
ACCESSIONS	BC113407
VECTOR	pDONR223
NCBI SYMBOL	UACA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]