



Reference: PTXBC113443

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SAPS1; KIAA1115; PP6R1; SAP190; serine/threonine-protein phosphatase 6 regulatory subunit 1; SAPS domain family, member 1; protein phosphatase 6 regulatory subunit 1		

For research use only.

GENE DESCRIPTION	SAPS domain family, member 1
GENE ID	22870
ORIGIN SPECIES	Human
ACCESSIONS	BC113443
VECTOR	pDONR223
NCBI SYMBOL	SAPS1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]