



Reference: PTXBC113454

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	U2AF1-RS2; U2AF1L2; U2AF1RS2; URP; ZC3H22; U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 2; U2 small nuclear ribonucleoprotein auxiliary factor, small subunit 2; U2(RNU2) small nuclear RNA auxiliary factor 1-like 2; U2		

GENE DESCRIPTION	zinc finger (CCCH type), RNA-binding motif and serine/arginine rich 2
GENE ID	8233
ORIGIN SPECIES	Human
ACCESSIONS	BC113454
VECTOR	pDONR223
NCBI SYMBOL	ZRSR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1