



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

HNRNPA3-heterogeneous nuclear ribonucleoprotein A3 Gene

Reference: PTXBC113470

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	2610510D13Rik; FBRNP; HNRPA3; heterogeneous nuclear ribonucleoprotein A3		
For research use only.			

GENE DESCRIPTION	heterogeneous nuclear ribonucleoprotein A3
GENE ID	220988
ORIGIN SPECIES	Human
ACCESSIONS	BC113470
VECTOR	pDONR223
NCBI SYMBOL	HNRNPA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]