

LOC100131551-hypothetical LOC100131551 Gene

Reference: PTXBC113537

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC100131551
GENE ID	100131551
ORIGIN SPECIES	Human
ACCESSIONS	BC113537
VECTOR	pDONR223
NCBI SYMBOL	LOC100131551
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaaagaatccaaggactgtgtcgtgggtgtcagctgctattctgtgccctggggacgaagaaggatctggagggaagactgtcacctcttgcgaatgtttccaagaagcagctcagcatgggagagaaaagggagatccagttccaggaacacctggccagtggttcacc
tggttatcgagtgacaggaccagccagcccttcccacctctcccacatgatgatggagcctggccttggagattcctgcgaaggaaatga