

PRY-PTPN13-like, Y-linked Gene

Reference: PTXBC113548

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRY1; PTPN13LY; PTPN13-like protein, Y-linked; testis-specific PTP-BL-related protein on Y; PTPN13-like, Y-linked		

For research use only.

GENE DESCRIPTION	PTPN13-like, Y-linked
GENE ID	9081
ORIGIN SPECIES	Human
ACCESSIONS	BC113548
VECTOR	pDONR223
NCBI SYMBOL	PRY
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagccactgggcttggcttctacttctctggagacaagacaatttgaatggcactgactgccaggatgcaatatttatacttctctgagactacggggagcatgtgttctgaacttcttgaacagaggtcttgaggccagaaggaagagatctaaagactcatttctctggagat
atgggaagggtggctgtatctcacttccactctgtgagatgaccgctggattaacccaccccaatttcagagatttccaaggctaccaccagagggtgcacggagctgatgcactgagcctgcaaaccaactctctgagaagcagggttatcttcacagtgctcggacagagcttctct
caggacactcgagagagccgtgttccaggcacttggggacatctgtggccactgtcatgaagaagactaa