



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C16orf72-chromosome 16 open reading frame 72 Gene

Reference: PTXBC113570

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0472 protein C16orf72; PRO0149; chromosome 16 open reading frame 72		
For research use only.			

GENE DESCRIPTION	chromosome 16 open reading frame 72
GENE ID	29035
ORIGIN SPECIES	Human
ACCESSIONS	BC113570
VECTOR	pDONR223
NCBI SYMBOL	C16orf72
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagagacgcggaagaggagagcgagcgccagatcaggagacagcgcccgagcactggtctccaagtgggagcgcagctcgtccgagcgcaaacgacgagcagctgcccccagctgcaggaggagcgcggcgccgcgcagccgacacaagcag
cagaagctgtggtcactcttcagaactcggccaccgcgtgggccagctctacaagaccgagtgtgtcagcagccaggacttctctgtgggtccccctccaaaacgagcaccgccgtcaccactctcaaaagaaagcgtggataccatcaacgaagtgttatgttgaaatca
gattgcctatcagcagcacaataagattgtgttgcgttgggttaaaaaacgcagaagaactatctgtcgagaagattgatcagcttctgtgtgaaagttctccaccacgaactctagatctccccagaactgactgtatgttccctaacccgactcttcacggaactgactca
ctgttagagactgatttgcgaacctccggcagcgcatagctctgtcgtttagtgggtgcaatggctatgataagcgtgcgttcagaccacccgactctctacacatgaagcagtgatcgaatgctagtcgaagggagaatggactcatgatgtcattgaacattcatatcaga
gaaatgtggtgactgactctgttgcactagaagacattctcagcccgatctggtgcgtatgctcagccacataaacgcacaagatgactcaa

More info: Please contact us for any question: contact@proteogenix.fr