



Reference: PTXBC113577

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	tolloid-like protein 2; tolloid like 2		
For research use only.			

GENE DESCRIPTION	tolloid-like 2
GENE ID	7093
ORIGIN SPECIES	Human
ACCESSIONS	BC113577
VECTOR	pDONR223
NCBI SYMBOL	TLL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]