



Reference: PTXBC113609

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CNR1; CNRN1; CRNR1; PCDH-ALPHA4; protocadherin alpha-4; KIAA0345-like 10; PCDH-alpha-4; ortholog of mouse CNR1; protocadherin alpha 4		

For research use only.

GENE DESCRIPTION	protocadherin alpha 4
GENE ID	56144
ORIGIN SPECIES	Human
ACCESSIONS	BC113609
VECTOR	pDONR223
NCBI SYMBOL	PCDHA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]