

PDCD1LG2-programmed cell death 1 ligand 2 Gene

Reference: PTXBC113678

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B7DC; Btdc; CD273; PD-L2; PDCD1L2; PDL2; bA574F11.2; programmed cell death 1 ligand 2; B7 dendritic cell molecule; B7-DC; PD-1-ligand 2; PDCD1 ligand 2; butyrophilin B7-DC; programmed death ligand 2		

For research use only.

GENE DESCRIPTION	programmed cell death 1 ligand 2
GENE ID	80380
ORIGIN SPECIES	Human
ACCESSIONS	BC113678
VECTOR	pDONR223
NCBI SYMBOL	PDCD1LG2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatctctcctgctaattgtgagcctggaattgcagcttcaccagatagcagctttattcacagtgcagctccctaaggaaactgtacataatagagcatggcagcaatgtgaccctggaatgcaactttgacactggaagtcagtgaaccttgagcaataacagccagttgcaaaaggt
ggaaaatgatacatccccacaccgtgaaagagccactttgtctggaggagcagctgcccctagggaaggcctctgtccacatactcaagtccaagtggaggacgaaggacagtaccaatgcataatcatctatgggtgcctgggactacaagtcactgactctgaaagtcaggcttc
ctacaggaaaataaacactcacatcctaaggttcagaaacagatgaggttagagctcactgcccaggctacaggtttatcctctggcagaagtatcctggccaacgtcagcgttctgccaacaccagccactccaggacccctgaaggcctctaccaggtcaccaggtgttctgcgccta
aagccacccctggcagaaacttcagctgtgtgttctggaatactcacgtgagggaaactttggccagcattgacctcaaaagtcagatggaaccaggacccatccaacttggctgcttcacatttccctcctgcatcattgtttcatttcatagccacagtgatagccctaagaaa
acaactctgcaaaagctgtattcttcaaaagacacaacaaaaagacctgtcaccacaacaaagagggaagtgaacagtgtctatctga