



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CD79A-CD79a molecule, immunoglobulin-associated alpha Gene

Reference: PTXBC113731

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD79a molecule; CD79a molecule, immunoglobulin-associated alpha; CD79a antigen (immunoglobulin-associated alpha); IGA; MB-1; B-cell antigen receptor complex-associated protein alpha chain; MB-1 membrane glycoprotein; ig-alpha; membrane-bound immunoglobuli		

For research use only.

GENE DESCRIPTION	CD79a molecule, immunoglobulin-associated alpha
GENE ID	973
ORIGIN SPECIES	Human
ACCESSIONS	BC113731
VECTOR	pDONR223
NCBI SYMBOL	CD79A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctcgtggggctcagagagctctccaaagctgctcgtccacatctctctctctctcgtcgtcgtcgtctacctgggcccctgggtgcaggccctgfgagtcacaaaggtccagcatcatgatggtgcagctgggggaagacgccacttccaatgccgcacaatagcagcaaaa
cgcaaacgtcactcgtggggcgcgtctctcatggcaactacacgtggccccctgaggtctctgggccccggcgaggaccaccaatggttagctgtgatcatcagaatgtgaacaagagccatgggggcatalatagctgtgcggggctcggagaggggaacgagtcataccagcagctctcgtgg
cactcatctccgtgctgcgcagccgccccccagccctctcctgacatgggggcagggggcaggacgaaccgaatcatcacacgcaggagggatcatcctgttctcgtggctgtgctgctggcgtgctgctgacagctcgtcgtctgacagacgctgctgtgacagacgctgta
gtaatacctgaagatgaataaacctttatgaagcctgaacgtgcacatgctcctcatgtagtgaagcaatctccgcctcgcagccactcaacagatgagctggcgaacctccaatcatalagagatggtccagctgaggaagccctga

More info: Please contact us for any question: contact@proteogenix.fr