



Reference: PTXBC113953

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HTRP-1; TRP1; short transient receptor potential channel 1; TRP-1; capacitative calcium channel protein Trp1; transient receptor potential canonical 1; transient receptor protein 1; transient receptor potential cation channel subfamily C member 1		

For research use only.

GENE DESCRIPTION	transient receptor potential cation channel, subfamily C, member 1
GENE ID	7220
ORIGIN SPECIES	Human
ACCESSIONS	BC113953
VECTOR	pDONR223
NCBI SYMBOL	TRPC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]