



Reference: PTXBC113978

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPR64; CBAVDX; EDDM6; TM7LN2; adhesion G-protein coupled receptor G2; G protein-coupled receptor 64; G protein-coupled receptor, epididymis-specific (seven transmembrane family); epididymal protein 6; human epididymis-specific protein 6; adhesion G protei		

GENE DESCRIPTION	G protein-coupled receptor 64
GENE ID	10149
ORIGIN SPECIES	Human
ACCESSIONS	BC113978
VECTOR	pDONR223
NCBI SYMBOL	GPR64
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]