

CTAG2-cancer/testis antigen 2 Gene

Reference: PTXBC113998

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAMEL; CT2; CT6.2; CT6.2a; CT6.2b; ESO2; LAGE-1; LAGE2B; cancer/testis antigen 2; CTL-recognized antigen on melanoma; LAGE-1a protein; autoimmunogenic cancer/testis antigen NY-ESO-2; cancer/testis antigen 6.2; cancer/testis antigen family 6, member 2a; ca		

For research use only.

GENE DESCRIPTION	cancer/testis antigen 2
GENE ID	30848
ORIGIN SPECIES	Human
ACCESSIONS	BC113998
VECTOR	pDONR223
NCBI SYMBOL	CTAG2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggccgaagccggggcacagggggttcgacggcgatgctgatggccaggagccctggcattcctgatggccaggggcaatgctggcgccaggagagggcgggtgccacggggcgagagggtccccggggcgaggggcagcaaggccctcggggccga
gaggaggcgccccgggggtccgcatggcggtgcccgttctgcgaggatggaaggtgccccctgcggggccaggagggccgacagccgctgcttgagttgcacatcacgatgcctttctcgcgccatggaagcgagctgtccgaggaatcctgccccgatccgcacc
gctccccgaccaggggcggttctgaaggacttcaccgtgtccggcaacctactgtttatccgactgactgtgcagaccaccgcaactgcagctctccatcagctcgtgtctccagcagctttccctgttgatgtggatcacgcagtgctttctgcccgtgttttggctcaggctccctcagg
gcagaggcgctaa