

RAB5C-RAB5C, member RAS oncogene family Gene

Reference: PTXBC114439

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB5C, member RAS oncogene family; RAB5C, member of RAS oncogene family; L1880; RAB5CL; RAB5L; RABL; ras-related protein Rab-5C		

For research use only.

GENE DESCRIPTION	RAB5C, member RAS oncogene family
GENE ID	5878
ORIGIN SPECIES	Human
ACCESSIONS	BC114439
VECTOR	pDONR223
NCBI SYMBOL	RAB5C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggtcggggaggcgcagcagacccaatggaccagctgctgggaacaagatctgtcaatttaagctggttctgctgggggagctctgcggtaggcaaatccagcctcgtcctccgtttgtcaaggagacagtttcacgagtaccaggagagacacaattggagcggccttcctcaca
cagactgtctgctggtgatgacacacagtcagtttgagatctgggacacagctggacagagcgggtatcagcctggccccatgtactatcggggggccaggtgccatcgtgtgtctatgacatcaccaacacagataattgcacgggccaagaactgggtgaaggagctaca
gaggcaggccagccccaacatcgtcattgcactcggggtacaaggcagacctggccagcaagagagccgtggaattccagggaagcacaagcctatgcagacgacaacagtttgctgttcatggagacatcagcaagactgcaatgaacgtgaacgaaatcttcagcaatagct
aagaagcttcccaagaacgagccccagaatgcaactggtgctcaggccgaaaccgaggtgtggacctcagagagaacaaccagcagcgagccagctgctgcagcaactga