



Reference: PTXBC114493

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	JNKBP-1; JNKBP1; mitogen-activated protein kinase-binding protein 1; JNK-binding protein 1; mitogen-activated protein kinase binding protein 1		

For research use only.

GENE DESCRIPTION	mitogen-activated protein kinase binding protein 1
GENE ID	23005
ORIGIN SPECIES	Human
ACCESSIONS	BC114493
VECTOR	pDONR223
NCBI SYMBOL	MAPKBP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr