



Reference: PTXBC114505

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GALNAC-T6; GalNAcT6; polypeptide N-acetylgalactosaminyltransferase 6; GalNAc transferase 6; UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 6; UDP-N-acetyl-alpha-D-galactosamine: polypeptide N-acetylgalactosaminyltransferase 6; UDP-N-acetyl-alpha		

For research use only.

GENE DESCRIPTION	UDP-N-acetyl-alpha-D-galactosamine:polypeptide acetylglactosaminyltransferase 6 (GalNAc-T6)	N-
GENE ID	11226	
ORIGIN SPECIES	Human	
ACCESSIONS	BC114505	
VECTOR	pDONR223	
NCBI SYMBOL	GALNT6	
LENGTH	ProteoGenix bp	
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG	

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr