

## KIF1B-kinesin family member 1B Gene

Reference: PTXBC115395

|                |                                                                                                                                   |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                               | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                              | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | kinesin superfamily protein KIF1B; kinesin-like protein KIF1B; CMT2; CMT2A; CMT2A1; HMSNII; KLP; NBLST1; kinesin family member 1B |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | kinesin family member 1B                                                                         |
| <b>GENE ID</b>          | 23095                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC115395                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | KIF1B                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgtcgggagccctcagtgaaagtggctgtccggtaaggccctcaattctcagagaccagcaaggaatccaaatgcattcagatgcaagcgacctcgaccagattattataccaaagaatccaaaggagctccaaagtccctcagcttcgacttctactgtctcattacac  
ccgaagatccctgtttgcatctcaaaacggtgtgacaatgacattggcaaggaaatgctctacacgctttgaggatataatgtctgtattttgcctatgggcagactggtgtgctggaatcttatacaatgatgggtaacaagaagaagccaggtgcatctccacagtattgtga  
agaactttttgagaaaatcaatgacaactgtaataagaatgtcttactgtgaggggtgagctacatggaaatftactgtgaagagtagcagagatttgcctgaatccaaaacaagggttaatttgcgtgtgctgtgaacacccactcttggacctatgtggagatctgccaagtggcagtt  
acttctacacagacalltgcctcagctcagtgatgctgggaacaaagccaggacagtgccagctacaaatgaatgaacaagtagccgttccacgctgtgtttacgattgttttaccagaagaacacagataatgagaccaacttccactgagaaggtcagtaaatcagcttggtg  
gatctagcaggagtgacagcagctgattcaactgggtccaaaggagctcattaaaggagagcaaatattataatgtcttacaactttggcgaagctattcagccttggccagggtgagtaaaagaagaagaacagattttatccctacagggtatgtacttactgtgctcctt  
cgagaaaatttaggtgcaattctcggactgcaatggtgtgctgtcagcccccggaatacaactcagatgagactttgagcactctgagatagcagatcgtgcaaaacaattaaatgcaatgctgttataatgaggaccccaatgccaaactgtgtcgtgaattaaaggagaggtga  
cacggctgaaggaccttctcgtgctcagggcctgggagatattattgatacatccatggggctccctcactcattcccatcttctgctcactcagtagtcagggtggcgttgacgtctgtgaccagatttcaagagagagatctgtacacctggagagaggaagctattgaacgtttaag  
gaatcagagaagatcattgctgagttgaatgaaactgggaagagaagcttgaataacagagggccatcagaatggagagagaggttggctgagatgggagtgccattcgggaagatgggaagaccttgggttttctacctaagaagaccccatctgttaacctcaatga  
agaccactaattgtctgagtgccacttattatcatcaaatggaatfacaagggttggccagcagatgctgagcgcgccaggacatagtgctgagcggggctcacattaagaagagcattgtatctccgagtgagagaagcagcggggaggtatctgtgaccttagagcccc  
tgtgagcgtcagaaaactacgttaattggcaagagggtgtccagcctgttcagctgcgctcaggaacacgtatcatcatgggtaaaaacattgtttccgctttaaccaccgggaacagcagcagctgagcagagagaagactcttctgctgagacccctctgagcctgtggaactgga  
catttgcagagggagcttctggaatacaaggaaatgatatgaacaagagatgggaagagggctacaggaatggagatcttatacaaaaaggagaaggaagacagatctcttttggagcagcagagactggacgcggtattctgatacgggggacgattctgacaagaggtgc  
tgtgaagagagctggaaactgattactctctgagagaaaagctacctccagcaagttgcaaacattgtttaaaaatgtgctcccaagcagtggaagaagacgtgaaccaatataatgtatcagatacccccaagaagggcgttgatgaagattcgaagtggtgcacaatcagat  
cttaaaattcaggctgcaaaagagatttgcctgaggttgcctcattcaggtcaggtcaggtcggcagagattgaagccctggccattgtcaagatgaaggagcttggcattgtatggcaagaagaccccaatgagcgggactctgagggcagtgccaggagctctgggata  
ccgtcgtgttggggatgagaagatcgaaagcgtcatggccactgggaaggccagcactgtatgtagactcaaggttcatatagacaagctggaagatattttgcaagaagtcaaaagcaaaataacatgaaagacagggagataaaagtcttaagaataaaatgctcaaaatgg  
aaaaagcttctgcaactgaggtcaggaacagaaaagccaggaagccacaaaacaaaggagcctgttggctgtgtgttagtagcaccttgagataatgtaagtaaaaggagacaatggagaactgcaaaagaagaacgtgtttccagctgtatgaatgggagtcagcttttag  
acgtggacgtctgcctgtgagtgagcgaagacaatttgcgtttaagaacttgcaacagcaggagatacaaaagcagcttctgctgcgcaagattgtacctataggttcatccctctgagaacccggaagccccgcttccctttaaagagcaacctaaacacagaacacttggagctctggg  
acacatacatcataacagaagatgaggtttatagagcttaggattccaaagacgatgaagcaggaaggaataagaagagagagcaagaaaaaggggttaaggagcttttaaggatccccagtttccatggggtctcaaggaatgagaagtcagatcacatccaagttagcaa  
gcagcaccattataatcagcaacagccactcactactgttggaagcaattctctcaataatggccagccgaaaagtagcgcgtccagagcatctgctcccgaggtagcttaaaactccacagtggtgcccccactgctgatgtacagacttccaggcgaagcccatctacac  
accgtcagcttactgttaattataacactggagtcagttagaggcgaatgcagccacttctatcagaagcagactgacaacacccagccactgtagccagttgtgtacacctccggtgagagacagttctcagcccaatctcaaaagtgtgtgcagaacacacagataaa

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)