



Reference: PTXBC117125

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATX1; D6S504E; SCA1; ataxin-1; alternative ataxin1; spinocerebellar ataxia type 1 protein; ataxin 1		

For research use only.

GENE DESCRIPTION	ataxin 1
GENE ID	6310
ORIGIN SPECIES	Human
ACCESSIONS	BC117125
VECTOR	pDONR223
NCBI SYMBOL	ATXN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]