



Reference: PTXBC117127

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	JAW1-related protein MRVI1; protein MRVI1; JAW1L; IP3R-associated cGMP kinase substrate; inositol 1,4,5-triphosphate receptor-associated cGMP kinase substrate; inositol 1,4,5-trisphosphate receptor-associated cGMP kinase substrate; murine retrovirus integ		

For research use only.

GENE DESCRIPTION	murine retrovirus integration site 1 homolog
GENE ID	10335
ORIGIN SPECIES	Human
ACCESSIONS	BC117127
VECTOR	pDONR223
NCBI SYMBOL	MRVI1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]