

SETD2-SET domain containing 2 Gene

Reference: PTXBC117162

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	histone-lysine N-methyltransferase SETD2; HBP231; HIF-1; HIP-1; HSPC069; HYPB; KMT3A; LLS; SET2; p231HBP; huntingtin interacting protein 1; huntingtin yeast partner B; huntingtin-interacting protein B; lysine N-methyltransferase 3A; SET domain containing		

For research use only.

GENE DESCRIPTION	SET domain containing 2
GENE ID	29072
ORIGIN SPECIES	Human
ACCESSIONS	BC117162
VECTOR	pDONR223
NCBI SYMBOL	SETD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr