

RAB41-RAB41, member RAS oncogene family Gene

Reference: PTXBC117239

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB41, member RAS oncogene family; RAB41, member RAS homolog family; ras-related protein Rab-41		

For research use only.

GENE DESCRIPTION	RAB41, member RAS oncogene family
GENE ID	347517
ORIGIN SPECIES	Human
ACCESSIONS	BC117239
VECTOR	pDONR223
NCBI SYMBOL	RAB41
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgcttgggtcagcagcagggcctggatggaggcggaggccttgggtctggaggctgccgaaagaacggaatacagctctgtgcaaatctaaactcttattcctgggagagcagagcgggaagacatccatcagccgttcattacaacagcttcggctgcgctgccagg
caactgttggaattgacttctgtctaagaccatgtacttggaggacaaatagttcagctgcagctatgggacacagctggccaggagcgtttcacagcctaattcctagctacattcgtgattcaactattgcagtggtgtctatgacattacaacatcaattctttaaggagacagataagt
gggtagaacacgctgcagcagaagagggtgacgatgtgtcatcatgttgggttaacaagattgattggataacaaaagacaagtcactgcagaacagggtgaagaaaaatccagaacctaattgtatgttattgagaccagtgcctaaacgggttacacgtgaaaaagctgttc
cggcgtgtgcttctgccctcttccacaaggacttcactccacaaaagaggggacgggtgaaatcgaactggaatccttcaggagtcaggcaacagaagctattgtga