



Reference: PTXBC117289

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	0610011B16Rik; CRN7; POD1; coronin-7; 70 kDa WD repeat tumor rejection antigen homolog; coronin 7		

For research use only.

GENE DESCRIPTION	coronin 7
GENE ID	79585
ORIGIN SPECIES	Human
ACCESSIONS	BC117289
VECTOR	pDONR223
NCBI SYMBOL	CORO7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1