



Reference: PTXBC117316

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	hGCMb; chorion-specific transcription factor GCMb; GCM motif protein 2; glial cells missing 2; glial cells missing homolog b; glide/gcm protein homolog; glial cells missing homolog 2		

For research use only.

GENE DESCRIPTION	glial cells missing homolog 2 (Drosophila)
GENE ID	9247
ORIGIN SPECIES	Human
ACCESSIONS	BC117316
VECTOR	pDONR223
NCBI SYMBOL	GCM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr