

CRYGD-crystallin, gamma D Gene

Reference: PTXBC117338

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CACA; CCA3; CCP; CRYG4; CTRCT4; PCC; cry-g-D; gamma-crystallin D; gamma crystallin 4; gamma-D-crystallin; crystallin gamma D		

For research use only.

GENE DESCRIPTION	crystallin, gamma D
GENE ID	1421
ORIGIN SPECIES	Human
ACCESSIONS	BC117338
VECTOR	pDONR223
NCBI SYMBOL	CRYGD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaagatcaccctctacgaggaccggggtccaggccgccactacgaatgcagcagcgaccacccaacctgcagccctacttgagccgctgaactcggcgcgctggacagcggctgctggtatgctctatgagcagcccaactactcgggctccagttctctgcgc
cgcgcgactatgccaccaccagcagtgatgggctcagcagctcgtccgctcctcggcctcaccctcctcgtcctcagcaggtacagactctatgagaggagggactacagagccagatgatagttactgaggactgctcctgtcttcaggaccgttcgcttcaatg
aaatccactccctcaacgtgctggaggctcctgggtcctctacgagctgtccaactaccgaggacggcagttacctgctgatgccaggggactataggcgctaccaggactggggggccacgaatgccagagtgggctctctgaggagagtcataattctcctga